




LAB RETINAL IMPLANTS

Labtician Ophthalmics **Retinal Implants** Catalogue




LAB RETINAL IMPLANTS

Lincoff rules for finding primary break




Superior Temporal or Nasal Detachments

In 98% of the cases, the primary break lies within 1-1/2 clock hours of the highest border.



Total or Superior Detachments that Cross the 12 o'clock Meridian

In 93% of the cases, the primary break is at 12 o'clock or in a triangle, the apex of which is at the ora serrata, and the sides of which intersect the equator 1 hour on either side of 12 o'clock.



Inferior Detachments

In 93% of the cases, the higher side of the detachment indicates on which side of the 6 o'clock radian an inferior break lies. When an inferior detachment is bullous, the primary break lies above the horizontal meridian.

Cross reference chart

Style	Product	Implant Type	Style	Product	Implant Type
32	S 2968	Grooved Strip	503	S 1982-3	Silicone Sponge Round
40	S 2969	Circling Band	504	S 1982-4	Silicone Sponge Round Silicone
➤ 40	S 5020	Silicone Lace	505	S 1982-5	Sponge Round
41	S 2970	Silicone Strip	506	S 1981-5	Silicone Sponge Oval
42	S 2971	Silicone Strip	506 G	S 1983-3	Silicone Sponge Grooved
70	S 3018	Silicone Sleeve	508 G	S 1983-4	Silicone Sponge Grooved
72	S 3071	Silicone Sleeve	509 G	S 1983-3.5	Silicone Sponge Grooved
219	S 2992	Grooved Strip	510	S 1984-5	Silicone Sponge Part Thickness
220	S 2998	Grooved Strip	➤ 510 G	S 1983-2.3	Silicone Sponge Grooved
225	S 3013	Grooved Strip	511	S 1984-7.5	Silicone Sponge Part Thickness
240	S 2987	Circling Band	515	S 1984-3.3	Silicone Sponge Part Thickness
➤ 240	S 5025	Silicone Lace	516 G	S 1983-2.3	Silicone Sponge Grooved
270	S 3019	Silicone Sleeve	907	S 1981-7	Silicone Sponge Oval
275	S 2994	Silicone Tire	➤ 2950	S 2950	Circling Band
276	S 2999	Silicone Tire	3083	S 3083	Silicone Sleeve Oval
277	S 2986	Silicone Tire	3084	S 3084	Silicone Sleeve Oval
279	S 2991	Silicone Tire	➤ 3093	S 3093	Silicone Sleeve Buckling
280	S 2996	Silicone Tire	4050	S 4050	Silicone Strip Wide
286	S 3010	Silicone Tire	➤ 8270	S 1982-7	Silicone Sponge Round
287	S 3014	Silicone Tire	➤ 8425	S 1984-2.5	Silicone Sponge Part Thickness
287 WG	S 3014 L	Silicone Tire	➤ 8570	S 1985-7	Silicone Sponge Scholda
501	S 1981-4	Silicone Sponge Oval			
502	S 1982-2	Silicone Sponge Round			

Labtician Ophthalmics, Inc.

2140 Winston Park Drive, Unit 6, Oakville, Ontario, Canada L6H 5V5

Telephone: 905.829.0055 Toll free: 1.800.265.8391 (Canada Only)

labtician.com

Order Desk & Customer Service

Toll free: 1.855.651.4934 (Canada & USA only)

Toll free fax: 1.855.390.9790 (Canada & USA only)

Email: orders@labtician.com

Manufactured by
LABTICIAN™
Moving Vision Forward

➤ Indicates products exclusive to Labtician

LAB1017 REV B

The **widest range** of scleral buckling implants available



CE 2797

Circling Bands

Style	Product
40	S 2969

To be used with 1, 2, 4

Style	Product
240	S 2987

To be used with 1, 3, 4

Style	Product
2950	S 2950

To be used with 1, 3, 4

Wong Buckling SleeveTM

Style	Product
3093	S 3093

4

Silicone Sleeves

Style	Product
70	S 3018

2

Style	Product
270	S 3019

3

Style	Product
72	S 3071

6

Style	Product
3083	S 3083

1

Style	Product
3084	S 3084

5

Silicone Strips

Style	Product
41	S 2970

To be used with 5, 6

Style	Product
42	S 2971

To be used with 5, 6

Silicone LaceTM

Style	Product
40	S 5020

To be used with 1, 2, 4

Style	Product
240	S 5025

To be used with 1, 3, 4

Grooved Strips

Style	Product
219	S 2992

To be used with 5, 6

Style	Product
32	S 2968

To be used with 5, 6

Style	Product
220	S 2998

To be used with 5, 6

Style	Product
225	S 3013

To be used with 5, 6

Wide Strips

Style	Product
4050	S 4050

To be used with 5, 6

Silicone Sponges Round

Style	Product	Style	Product	Style	Product	Style	Product	Style	Product
502	S 1982-2	503	S 1982-3	504	S 1982-4	505	S 1982-5	8270	S 1982-7

Silicone Sponges Grooved

Style	Product	Style	Product	Style	Product	Style	Product
506G	S 1983-3	508G	S 1983-4	509G	S 1983-3.5	510G	S 1983-2.5

Grooved

Silicone Sponges Part Thickness

Style	Product	Style	Product	Style	Product	Style	Product	Style	Product
516G	S 1983-2.3	8425	S 1984-2.5	515	S 1984-3.3	510	S 1984-5	511	S 1984-7.5

Silicone Sponges Oval

Style	Product	Style	Product	Style	Product
501	S 1981-4	506	S 1981-5	907	S 1981-7

Silicone Sponges Scholda

Silicone Tires Convex

Style	Product	Style	Product	Style	Product	Style	Product
8570	S 1985-7	286	S 3010	287	S 3014	287WG	S 3014L

Silicone Tires Concave

Silicone Tires Asymmetrical

Style	Product	Style	Product	Style	Product	Style	Product	Style	Product
277	S 2986	279	S 2991	275	S 2994	276	S 2999	280	S 2996